# Mastercycler pro Self Test Function

Instructions for Use  
Gebrauchsanweisung  
Notice d'utilisation  
Manual de instrucciones  
Manual d'uso

**eppendorf**

Copyright© 2008 Eppendorf AG, Hamburg. No part of this publication may be reproduced without the prior permission of the copyright owner.

**Trademarks**

eppendorf, Mastercycler and vapoprotect are registered trademarks of Eppendorf AG, Hamburg, Germany.

Registered trademarks are not marked in all cases with <sup>TM</sup> or <sup>®</sup> in this manual.

**Instructions for Use. . . . . 2**

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**Manual de instrucciones. . . . . 27**

**Manual d'uso . . . . . 35**

<b>1</b>	<b>User instructions</b>	<b>4</b>
1.1	Using this manual	4
<b>2</b>	<b>Product description</b>	<b>4</b>
2.1	Delivery package	4
2.2	Features	4
<b>3</b>	<b>Operation</b>	<b>4</b>
3.1	Activate self test function	4
3.2	Perform self test	5
3.2.1	Create a certificate	7
3.3	Display, print or export history	8
3.3.1	Print History	9
3.3.2	Export history as file	9
<b>4</b>	<b>Troubleshooting</b>	<b>9</b>
<b>5</b>	<b>Transport, storage and disposal</b>	<b>10</b>
5.1	Disposal	10
	<b>Ordering information</b>	<b>43</b>

# Mastercycler pro Self Test Function – Instructions for Use

EN

Instructions for Use

## 1 User instructions

### 1.1 Using this manual

These operating instructions for the self test function of the Mastercycler pro is a supplement to the Mastercycler pro operating manual and is no replacement.

You should therefore also read the operating manual before using the centrifuge for the first time. The current version of the instructions for use is available on the internet at [www.eppendorf.com](http://www.eppendorf.com)

## 2 Product description

### 2.1 Delivery package

Quantity	Order No. (International)	Order No. (North America)	Description
1	6320 071.001	950030040	Selftest USB Key for Mastercycler pro
1	6320 071.168	-	Mastercycler pro Selftest - Instructions for Use

### 2.2 Features

The self test function can be used with the Control Panel with USB connections combined with the following devices:

- Mastercycler pro
- Mastercycler pro S
- Mastercycler pro 384

The self test function permits a quick check regarding the following attributes of the Mastercycler pro thermoblock without requiring an external measurement system:

- Heating and cooling rates
- Correct function of the temperature control circuit
- Temperature homogeneity within each of the three block zones (left, middle and right zone)
- Temperature homogeneity throughout the whole block

After a passed self test a certificate in PDF format can be generated.

The self test USB key activates the self test function on the Control Panel for all connected Mastercycler pro.

## 3 Operation


### 3.1 Activate self test function

#### Requirement

- The Mastercycler pro is switched on.
  - A user is logged on (if the function **PIN** is activated).
  - The navigation tree is displayed on the Control Panel.
1. Connect the self test USB key to a free USB connection on the right side of the Control Panel.

## Mastercycler pro Self Test Function – Instructions for Use

The message **Please wait** appears. After a couple of seconds the self test USB key is identified by the Control Panel.

- Open the node of the Mastercycler pro in the navigation tree.

The self test function is now unlocked.

The status of the self test function can be identified by its symbols in the navigation tree.

- Self test function activated
- Self test function locked

The self test function is activated until the self test USB key is disconnected from the Control Panel.

### 3.2 Perform self test


#### **Warning! Risk of burns from the hot thermoblock and heated lid.**

The thermoblock and the heated lid very quickly reach temperatures above 95 °C.

- Keep heated lid closed until the temperature of the thermoblock is below 30 °C.

#### Requirement

- The Mastercycler pro is switched on and is in idle time.
  - A user is logged on (if the function **PIN** is activated).
  - The navigation tree is displayed on the Control Panel.
  - The self test function is unlocked (see p. 4).
- Equip the block of the Mastercycler pro with an empty PCR plate and close the heated lid.

- Mark the node **Self Test**.
- Press function key **Start**.

## Mastercycler pro Self Test Function – Instructions for Use

EN

Instructions for Use

Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	UNTESTED	UNTESTED	UNTESTED
Temperature control	UNTESTED	UNTESTED	UNTESTED
Temperature homogeneity	UNTESTED	UNTESTED	UNTESTED
Overall temp. homogeneity	UNTESTED		
Press "Start" to begin		Lid: --- °C	Block: --- °C
Start		Exit	

- Press function key **Start** again.

The self test runs automatically. Information on the currently running test is displayed.

Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	Test in progress...		
Temperature homogeneity			
Overall temp. homogeneity			
Determining block temperature homogeneity		Lid 105 °C Block 94.5 °C	
		Abort	

<b>Block zone</b>	Left, middle and right block zone according to the temperature control circuit of the block.
<b>Heating/cooling rates</b>	Temperature control velocity of the block during heating or cooling.
<b>Temperature control</b>	Accuracy of the thermoregulation.
<b>Temperature homogeneity</b>	Temperature homogeneity within each individual control circuit.
<b>Overall temp. homogeneity</b>	Temperature homogeneity across the entire block.

# Mastercycler pro Self Test Function – Instructions for Use

EN

Instructions for Use

After 5 to 10 minutes the test result is displayed.

Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	PASS	PASS	PASS
Temperature homogeneity	PASS	PASS	PASS
Overall temp. homogeneity	PASS		
TEST PASSED		Lid 105 °C Block 4.0 °C	
		Start Certificate Exit	

Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	PASS	FAIL	FAIL
Temperature homogeneity	PASS	FAIL	FAIL
Overall temp. homogeneity	FAIL		
TEST FAILED		Lid 105 C Block 3.9 C	
		Start Exit	

- Passed tests are labeled with **PASS** and highlighted in green.
- Failed tests are labeled with **FAIL** and highlighted in red.

The self test is passed if all subtests have been passed.

- i** Failure of a test points to alteration, maladjustment or malfunctioning of the block. In this case please contact Eppendorf Service.

- i** To verify and calibrate the temperature accuracy of the block use the Eppendorf temperature validation system (see on page 43).

## 3.2.1 Create a certificate

After the self test has been passed you can either export a certificate in PDF format onto a USB storage medium or export an MMC.

### Requirement

The self test has passed all subtests and the test result is displayed.

1. USB storage medium or MMC on the right side of the Control Panel.
2. Press function key **Certificate**

The certificate is copied as a PDF file onto the storage medium. After completing the procedure the path of the copied files is displayed.

You can display the certificate on a PC and print it. For this purpose you need the software Adobe Reader or another PDF display program.

The certificate includes fields for the validity date as well as the signature.

# Mastercycler pro Self Test Function – Instructions for Use

EN

Instructions for Use


Fig. 1: Example of a self test certificate.

## 3.3 Display, print or export history

The application of the self test function and the result are recorded for each Mastercycler pro in a log file. In order to display, print or export these to a storage medium proceed as follows.

### Requirement

- A user is logged on (if the function **PIN** is activated).
  - The navigation tree is displayed on the Control Panel.
- Open the node of the Mastercycler pro and mark the entry **Self Test**.
  - Press the function key **History**.

The self test history is displayed.

## Mastercycler pro Self Test Function – Instructions for Use

### 3.3.1 Print History

#### Requirement

A USB supported printer or a parallel printer is connected (see Mastercycler pro operating manual).

- ▶ Press function key **Print**.

### 3.3.2 Export history as file

1. Connect USB storage medium or MMC to the Control Panel.
2. Press function key **Export**.
3. Select file format PDF or TXT and confirm with **Ok**.

The file is copied to the storage medium. After completing the procedure the path of the copied files is displayed.

## 4 Troubleshooting

Symptom / message	Cause	Remedy
self test USB key is not identified by Control Panel.	The self test USB key was connected to an unsupported USB hub.	▶ Connect the self test USB key directly to the Control Panel.
	The USB connections of the Control Panel are not active because the Control Panel was connected to a Mastercycler ep (the symbol  is displayed).	▶ Connect the Control Panel directly to a Mastercycler pro (see Mastercycler pro operating manual)
	The USB connections were switched off due to a damaged or incompatible USB device (the symbol  is displayed).	▶ Disconnect all USB devices from Control Panel and restart them.
The self test node is not displayed	No user is logged on.	▶ Log on as user or as administrator.
	A device of the Mastercycler ep series is selected.	<ul style="list-style-type: none"> <li>▶ Select a Mastercycler pro node.</li> <li>▶ Mini satellite system: check the address setting of all connected cyclers (see Mastercycler pro operating manual).</li> </ul>

### 5 Transport, storage and disposal

#### 5.1 Disposal

In the event of disposing of the product, please observe the applicable legal regulations.

**Information on the disposal of electrical and electronic devices in the European Community:**

The disposal of electrical devices is regulated within the European Community by national regulations based on EU Directive 2002/96/EC pertaining to waste electrical and electronic equipment (WEEE).

In accordance with this, any devices delivered after 13/08/2005 on a business-to-business basis, which includes this product, may no longer be disposed of in household waste. To document this they have been marked with the following identification:

Because disposal regulations may differ from one country to another within the EU please contact your supplier if necessary.

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Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	UNTESTED	UNTESTED	UNTESTED
Temperature control	UNTESTED	UNTESTED	UNTESTED
Temperature homogeneity	UNTESTED	UNTESTED	UNTESTED
Overall temp. homogeneity	UNTESTED		
Press "Start" to begin		Lid: --- °C	Block: -- °C
		Start	Exit

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Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	Test in progress...		
Temperature homogeneity			
Overall temp. homogeneity			
Determining block temperature homogeneity.	Lid 105 °C Block 94.5 °C		
		Abort	

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DE

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Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	PASS	PASS	PASS
Temperature homogeneity	PASS	PASS	PASS
Overall temp. homogeneity	PASS		
TEST PASSED		Lid 105 °C Block 4.0 °C	
		Start	Certificate
			Exit

Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	PASS	FAIL	FAIL
Temperature homogeneity	PASS	FAIL	FAIL
Overall temp. homogeneity	FAIL		
TEST FAILED		Lid 105 C Block 3.9 C	
		Start	Exit

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Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	UNTESTED	UNTESTED	UNTESTED
Temperature control	UNTESTED	UNTESTED	UNTESTED
Temperature homogeneity	UNTESTED	UNTESTED	UNTESTED
Overall temp. homogeneity	UNTESTED		
Press "Start" to begin		Lid: --- °C	Block: --- °C
		Start	Exit

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Cycler1 Cycler Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	Test in progress...		
Temperature homogeneity			
Overall temp. homogeneity			
Determining block temperature homogeneity	Lid 105 °C Block 94.5 °C		
		Abort	

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Cyc1er1 Cyc1er Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	UNTESTED	UNTESTED	UNTESTED
Temperature control	UNTESTED	UNTESTED	UNTESTED
Temperature homogeneity	UNTESTED	UNTESTED	UNTESTED
Overall temp. homogeneity	UNTESTED		
Press "Start" to begin		Lid: --- °C	Block: --- °C
		Start	Exit

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Cyc1er1 Cyc1er Self Test			
Test Results	Block zone		
	Left	Center	Right
Heating/cooling rates	PASS	PASS	PASS
Temperature control	Test in progress...		
Temperature homogeneity			
Overall temp. homogeneity			
Determining block temperature homogeneity	Lid 105 C Block 94.5 C		
		Abort	

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Vaudaux-Eppendorf AG
Eppendorf UK Limited
Eppendorf North America, Inc.
Eppendorf AG · 22331 Hamburg · Germany · Tel. +49 40 5 38 01-0 · Fax +49 40 5 38 01-556
Eppendorf North America, Inc. · One Cantiague Road, P.O. Box 1019 · Westbury, N.Y. 11590-0207 USA